



LifePort™

Kidney Transporter

Perfusion Laboratory Recommendations

Lab Space (12' x 12' minimum)

Countertop with overhead cabinetry or shelving

Laminar Flow Hood (6')

Thermo Forma model #1284

Scrub Sink (25" x 22" x 10" or larger) with Eyewash Station

Kohler Ballad single-basin utility sink or Franke 20-gauge double-basin Delta model #54T5437A-R8 gooseneck faucet with double foot control

Backtable (18"-24" W x 4'-6' L)

Adjustable height, stainless

Stool

Wheeled, adjustable height

Ice Machine (50-70 lbs/day)

Crushed or pelleted ice preferred, near sink

Refrigerator/Freezer (10-18 cu ft)**Storage Shelving (74" H x 48" W x 18"-24" D)**

Wire, stainless, solid surface for pump storage, locking casters

Miscellaneous Supplies

Nonsterile and sterile surgical gloves (Sm, Med, Lg, XL)

Surgical scrub caps, masks, gowns

Surgical scrub brushes

Eye protection

Sterile table covers and hood utility drapes

Sterile towels or gauze

O-silk ties

Disposable scalpels and biopsy needles (optional)

Sterile bag decanters

Sterile normal saline liter bags (for sterile ice)

Plastic, stackable storage bins

Sharps, biohazard, and standard waste containers and trash bags

Cleaning supplies

Services

Laundry (or access to scrubs)

Biohazard waste pickup and disposal